

# A Sample of TIMS DICOM System Sites in Europe



## **Royal Free Hospital - London, England, UK**

- Philips cardiac cath lab .....Convert to DICOM, send to PACS

## **Maelor Hospital - Wrexham, Wales, UK**

- Siemens Siregraph CF fluoroscopy .....Convert to DICOM, send to Insignia PACS
- Siemens Siregraph TOP fluoroscopy .....Convert to DICOM, send to Insignia PACS
- Siemens Polystar fluoroscopy .....Convert to DICOM, send to Insignia PACS
- Toshiba nuclear medicine .....Convert to DICOM, send to Insignia PACS
- Siemens mobile C-arm .....Convert to DICOM, send to Insignia PACS
- Philips mobile C-arm .....Convert to DICOM, send to Insignia PACS
- Thermal paper scanning .....Convert thermal paper scans to DICOM, send to Insignia PACS

## **Bronglais Hospital - Aberystwyth, Wales, UK**

- Mobile MR trailer .....Import studies via CD, send to Insignia PACS
- Toshiba portable ultrasound .....Convert to DICOM at remote sites, send to Insignia PACS when return to the hospital
- Toshiba CT .....TIMS Consultant for transmitting the live patient study to a radiologist's office for study input
- General study import.....Import studies from CD, send to Insignia PACS

## **Morriston Hospital - Swansea, Wales, UK**

- Philips Integris bi-plane cardiac cath lab .....Convert to DICOM, send to GE Centricity PACS

## **University Hospital - Cardiff, Wales, UK**

- Philips cardiac cath lab .....Convert to DICOM, send to GE Centricity PACS

## **Prince Philip Hospital - Carmarthenshire, Wales, UK**

- Ultrasound .....Convert to DICOM, send to GE Centricity PACS

## **Instituto Europeo di Oncologia - Milan, Italy**

- X-ray angiography .....Convert to DICOM, send to PACS

## **National General Hospital of Piraeus - Piraeus, Greece**

- Philips cardiac cath lab .....Convert to DICOM, record to CD/DVD

## **University Hospital Groningen - Groningen, Netherlands**

- Fluoroscopy .....Convert to DICOM, send to Rogan PACS

## **Bakulev Science Center of Heart & Vessel Surgery - Moscow, Russia**

- HP Sonos ultrasound .....Convert to DICOM, review in place, send to DICOM review station